

## CENTURY EVER-ON MARKER



### Directs of Use:

- Remove Cap. For cold start press the rubber ball/valve 4-5 times at the bottom for ink flow at the tip.
- Mark vertically at 90-degree angle & ensure the ball presses the surface for smooth ink flow.
- As pressure reduces, press the rubber valve 2-3 times to continue writing.
- Do not press the rubber bulb when not in use as it will lead to increase ink pressure.
- Clean the tip & put the cap back on when after each use.

### Advantages

- Century's Everon Marker is a unique marking device used in textile mills & factories for permanent marking on gray cloth & fabrics.
  - It has a rigid aluminum body making it difficult to break, puncture or leak.
  - It also has a rubber valve at the bottom for smooth ink flow. When the valve is pumped it compresses ink towards a fine ball point for smooth flow and easy writing
  - The viscous ink makes permanent, bright & opaque markings on all kinds of fabrics like cotton, denim, polyester, wool, silk, linen etc. and remains visible even after the dyeing & other finishing processes.
  - It is widely used in dyeing houses & textile mills in processes like dyeing, scouring, bleaching, hosiery, knitting & weaving.
  - Markings are acid alkali resistant, heat resistant & remain permanent even after dyeing & finishing processes.
- Everon Marker is designed whereby it can be used up to last drop of ink when compared to traditional tube type markers where almost 20% of ink goes waste or we are unable to use after squeezing, as it is difficult to even hold the tube It is also known as 'pump-type marker' or 'TEXPEN'.

### Caution:

- Do not fold the tube back & forth to avoid leakage and breakage.
- Harmful if swallowed or intentional inhalation.
- Avoid contact with eyes. Keep away from open flame.
- Please test before use.

### Packaging Information

Century's Ever-On Marker	Century's Ever-On Marker	60 Ml Ink Content	—
	Century's Ever-On Marker 1 pc	87 Grams	—
	Century's Ever-On Marker 15 pc pack (Plastic Box)	1.50 Kg	18.5 X 13 X 8.5
	Century's Ever-On Marker 180 pc carton (15 x 12 plastic box)	21 Kg	43 X 40.5 X 19

### Physical Properties:

Appearance	Thick coloured paste
Odour	Characteristic strong odour
Flash Point (°C)	-4.5°C
Boiling Range (°C)	77 - 112
Specific Gravity	1.00 ± 0.05
VOC (% by weight)	Max. 58
Drying Time (full dry)	Max. 5 Minutes

### Tip Sizes Available

- **Point No.:** 2 Mm Markings, 3 Mm Markings, 4 Mm Markings

### Colors Available & Product Code :

**YELLOW** **RED** **BLACK** **BLUE** **ORANGE** **GREEN** **WHITE** **PINK** **VIOLET**

- **Note:** Colors pink & violet is make to order

