



CENTURY TEX MARKER

Directs of Use:

- Open cap. Insert the aluminum key provided with the marker at the bottom end of the tube & start rolling upwards for pressure.
- Mark on 90 degree angle for smooth ink flow. Keep rolling the key as pressure reduces.
- Do not fold the tube back & forth as it could lead to puncture.
- Do not use on other surface other than fabrics.
- Put the cap back on after use.

Features

- Viscous ink makes Permanent, Bright & Opaque Markings on all type of fabric surfaces.
- Collapsible Aluminum body makes it easy to squeeze out optimum quantity of ink with the help of aluminum key provided with marker.
- Resistance to Acids, Alkalis & other chemical aggressions in finishing processes.
- Markings can withstand all textile sorting / sizing and finishing processes.
- Quick identification of lots upon fabrics dyed in dark shades also.
- Most economical method of marking on gray fabrics before processing.

Storage

- Store locked up. Keep away from heat, sparks and open flame. Store it in well ventilated, cool & dry place away from direct sunlight. Store it in original tightly closed container.

Disposal

- Collect and reclaim or dispose in sealed containers at licensed waste disposal site. Dispose of contents/ container in accordance with local / regional / national / international regulations.

Tip Sizes Available

- **Point No. 2:** Writes 2 mm wide for general purpose
- **Point No. 3:** Writes 3 mm wide, specially for knitted / hosiery / denim Fabric and Canvas type fabric
- **Point No. 4:** For raised fabrics

Caution:

- Do not fold the tube back & forth to avoid leakage and breakage.
- Harmful if swallowed or intentional inhalation.
- Avoid contact with eyes. Keep away from open flame.
- Please test before use.

Packaging Information

Century's Tex Marker	Century's Tex Marker	60 Ml Ink Content	—
Century's Tex Marker	Century's Tex Marker	75 Grams	—
Century's Tex Marker With Monocarton (1pc)	Century's Tex Marker With Monocarton (1pc)	86 Grams	17.6 X 3.6 X 2.6
Century's Tex Marker 10 Pc With Monocarton (Plastic Box)	Century's Tex Marker 10 Pc With Monocarton (Plastic Box)	1 Kg	18.5 X 13 X 8.5
Century's Tex Marker 15 Pc Packing Without Monocarton (Plastic Box)	Century's Tex Marker 15 Pc Packing Without Monocarton (Plastic Box)	1.25 Kg	18.5 X 13 X 8.5
Century's Tex Marker 160 Pc Carton (10 Pc X 16 Plastic Box)	Century's Tex Marker 160 Pc Carton (10 Pc X 16 Plastic Box)	18.5 Kg	57 X 41 X 19
Century's Tex Marker 180 Pc Carton (15 Pc X 12 Plastic Box)	Century's Tex Marker 180 Pc Carton (15 Pc X 12 Plastic Box)	15.550 Kg	43 X 41 X 19
Century's Tex Marker 240 Pc Carton (10 Pc X 24 Plastic Box)	Century's Tex Marker 240 Pc Carton (10 Pc X 24 Plastic Box)	25.65 Kg	44.5 X 39.5 X 38

Physical Properties:

Appearance	Thick coloured paste
Odour	Characteristic strong odour
Flash Point (°C)	-4.5°C
Boiling Range (°C)	77 - 112
Specific Gravity	1.00 ± 0.05
VOC (% by weight)	Max. 58
Drying Time (full dry)	Max. 5 Minutes

Colors Available:



Note: Colors violet is make to order

